

AD-A171 946

EXPERIMENTAL RESEARCH ON OPTOGALVANIC EFFECTS(U)  
WISCONSIN UNIV-MADISON DEPT OF PHYSICS J E LAWLER  
05 JUN 86 AFOSR-RR-86-0518 AFOSR-84-2298


1/1

UNCLASSIFIED

F/G 14/2

NL





Resolution Test Chart  
 100% Contrast, 100% Modulation Transfer Function

AD-A171 946

REPORT DOCUMENTATION PAGE		READ INSTRUCTIONS BEFORE COMPLETING FORM
1. REPORT NUMBER <b>AFOSR-TR- 86-0518</b>	2. GOVT ACCESSION NO.	3. RECIPIENT'S CATALOG NUMBER
4. TITLE (and Subtitle) Final Scientific Report for DOD Equipment Grant AFOSR 84-0298 under the URIP		5. TYPE OF REPORT & PERIOD COVERED 08/01/84 through 07/31/85
		6. PERFORMING ORG. REPORT NUMBER
7. AUTHOR(s) J. E. Lawler		8. CONTRACT OR GRANT NUMBER(s) AFOSR-84-0298 <del>AFOSR-84-0328</del>
9. PERFORMING ORGANIZATION NAME AND ADDRESS Department of Physics, University of Wisconsin 1150 University Avenue Madison, WI 53706		10. PROGRAM ELEMENT, PROJECT, TASK AREA & WORK UNIT NUMBERS W110DF, D3C1, A7
11. CONTROLLING OFFICE NAME AND ADDRESS Air Force Office of Scientific Research, NP Bolling Air Force Base, D. C. 20332		12. REPORT DATE June 5, 1986
		13. NUMBER OF PAGES 4
14. MONITORING AGENCY NAME & ADDRESS (if different from Controlling Office) AFOSR/NP, BICIG 410 Bolling AB D. C. 20332-6148		15. SECURITY CLASS. (of this report) Unclassified
16. DISTRIBUTION STATEMENT (of this Report)  Approved for public release, distribution unlimited		15a. DECLASSIFICATION/DOWNGRADING SCHEDULE
17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, if different from Report)  Approved for public release, distribution unlimited		
18. SUPPLEMENTARY NOTES		
19. KEY WORDS (Continue on reverse side if necessary and identify by block number)		
20. ABSTRACT (Continue on reverse side if necessary and identify by block number) The equipment acquired under a DOD University Research Instrumentation Program Grant is described. Applications of the equipment in ongoing DOD sponsored research on optogalvanic effects is reported.		

MIC FILE COPY

DD FORM 1 JAN 73 1473

EDITION OF 1 NOV 65 IS OBSOLETE

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

86 9 15 050

AFOSR-TR. 86-0518

Final Technical Report

to the

Air Force Office of Scientific Research

Grant AFOSR-84-0298

under the

University Research Instrumentation Program

Principal Investigator: J. E. Lawler

Approved for public release,  
distribution unlimited

AIR FORCE OFFICE OF SCIENTIFIC RESEARCH (AFSC)  
NOTICE OF TRANSMITTAL TO DTIC  
This technical report has been reviewed and is  
approved for public release IAW AFR 190-12.  
Distribution is unlimited.  
MATTHEW J. KERPER  
Chief, Technical Information Division

86 9 15 07

A-1

pattern and determine the laser wavelength. This wavemeter is well suited for use with pulsed or c.w. dye lasers. Other wavemeters based on Michelson interferometers which have moving mirrors and which involve sequential detection of fringes are not suited for use with pulsed lasers. Pulsed dye lasers typically have very low ( $\sim 10^{-7}$ ) duty cycles. We built a cart for the Fizeau wavemeter and associated computer and oscilloscope. It is now shared among three laser laboratories and is in steady use.

The second item we purchased was a Spex 1403 Double Spectrometer. A double spectrometer with large holographic gratings, such as this one, produces very good rejection of stray light (approximately  $10^{-14}$  when  $20 \text{ cm}^{-1}$  off the bandpass). This exceedingly high rejection is essential in a variety of fluorescence and scattering experiments in discharge plasmas. It makes it possible to see very weak signals against bright background light at nearby wavelengths. We have no other instrument which achieves comparable performance. The Spex 1403 Double Spectrometer has also been built into its own cart so that it can be shared among several laser laboratories.

The final item we purchased under this grant was a Spellman High Voltage Regulated Power Supply. This is a 3 kV supply which can deliver 0.5 Amp. It is regulated to 0.01%. We need it to run large stable discharges. Optogalvanic diagnostics achieve maximum sensitivity in stable well behaved discharges. The high degree of regulation is essential if we are to get the maximum amount of information from optogalvanic diagnostics.

In summary three pieces of capital equipment were purchased under this grant including: a Fizeau laser wavemeter, a double spectrometer, and highly regulated high voltage power supply. All are being used in ongoing research on optogalvanic effects.

References

1. D. K. Doughty and J. E. Lawler, Appl. Physics Lett. 42, 234 (1983).
2. D. K. Doughty and J. E. Lawler, Appl. Physics Lett. 45, 611 (1984).
3. D. K. Doughty, S. Salih, and J. E. Lawler, Physics Lett. 103A, 41 (1984).

END

DTIC

10-86